



$^{188}\text{W}/^{188}\text{Re}$ Generator

Production of high-energy beta emitting rhenium-188

Produced by OncoBeta® GmbH

OncoBeta® GmbH is a medical device and radiochemical company, specializing in the development and commercialization of state-of-the-art, innovative radioisotope therapies and products utilizing rhenium-188. OncoBeta® offers an innovative Skin Cancer Therapy targeting non-melanoma skin cancers and produces tungsten (wolfram)-188/rhenium-188 ($^{188}\text{W}/^{188}\text{Re}$) Generators for commercial use.

Schleißheimer Str. 91 | 85748 Garching n. Munich (Germany) | Tel: +49 (89) 32667330 | info@oncobeta.com

¹⁸⁸W/¹⁸⁸Re Generator

Not for Human Use.

The OncoBeta® ¹⁸⁸W/¹⁸⁸Re Generator is based on an alumina column. The longer-lived mother radionuclide ¹⁸⁸W is adsorbed on the column material while the daughter radionuclide ¹⁸⁸Re is selectively eluted from the generator with sterile isotonic solution as sodium perrhenate (Na¹⁸⁸ReO₄).

USAGE

The ¹⁸⁸W/¹⁸⁸Re Generator can be used for [radiolabelling reactions](#) or directly as a [high dose liquid radioactive source](#).

MAIN FEATURES

- › Available in [different sizes](#) upon request.
Up to 150 GBq (4 Ci) ¹⁸⁸W
- › High activity concentration.
Up to 9 GBq/ml ¹⁸⁸Re

SPECIFICATIONS

Column material	Aluminium oxide (Al ₂ O ₃) packed in quartz glass column
Eluent	Sterile 0.9 % saline solution
Chemical form	Sodium perrhenate (Na ¹⁸⁸ ReO ₄)
Appearance	Clear and colorless solution
Elution volume	5 - 14 ml (with small volumes not completed elution)
Elution yield	75 - 85 %
Radiochemical purity	> 98 % Re(VII)
Tungsten breakthrough	¹⁸⁸ W < 0.01 % in relation to ¹⁸⁸ Re at calibration time
Total gamma impurities	< 0.01 % in relation to ¹⁸⁸ Re at calibration time
Shelf-life	10 months
Specific activity	Product is carrier free
Generator size (nominal ¹⁸⁸W activity)	3.7- 150 GBq
Generator dimensions	Diameter 145 mm, height 210 mm
Generator weight	16 kg
Interface	Standard: Stainless steel needles for elution with evacuated vials. On request: Tubings with Luer-Lock connectors (inlet female, outlet male)

**This product is for laboratory research purposes only*

$^{188}\text{W}/^{188}\text{Re}$ Generator

Not for Human Use.


ELUTION

The **elution yield** of rhenium-188 is **75-85 %**, and it remains stable throughout the shelf-life of the generator.

Set up

- › The $^{188}\text{W}/^{188}\text{Re}$ generator must be eluted with sterile 0,9 % saline solution *only*.
- › The maximum time for the $^{188}\text{W}/^{188}\text{Re}$ generator to be left uneluted is three (3) weeks.

Typical $^{188}\text{W}/^{188}\text{Re}$ generator elution profile



Elution:
Vials

The inlet/outlet cover of the OncoBeta® $^{188}\text{W}/^{188}\text{Re}$ Generator offers two holes. The smaller IN hole with two needles for **saline vial** and the larger OUT hole with one needle for **evacuated vial**.

In the elution process, vacuum moves the saline solution through the column into the evacuated vial.

Elution:
Syringe, pump, dispenser system

On request, the OncoBeta® $^{188}\text{W}/^{188}\text{Re}$ Generator is delivered with tubing attached on the inlet and outlet needles that have **Luer Lock connectors**. By these means, the generator can be connected to a syringe, a peristaltic pump or a dispenser.



Produced by OncoBeta® GmbH

OncoBeta® GmbH is a medical device and radiochemical company, specializing in the development and commercialization of state-of-the-art, innovative radioisotope therapies and products utilizing rhenium-188. OncoBeta® offers an innovative Skin Cancer Therapy targeting non-melanoma skin cancers and produces tungsten (wolfram)-188/rhenium-188 ($^{188}\text{W}/^{188}\text{Re}$) Generators for commercial use.

Schleißheimer Str. 91 | 85748 Garching n. Munich (Germany) | Tel: +49 (89) 32667330 | info@oncobeta.com

Version C | Date: 22.02.2021